

DNAJC5B antibody

Product Information

Catalog No.:	FNab02467
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 5 beta
Synonyms:	DnaJ homolog subfamily C member 5B Cysteine-string protein isoform beta (CSP-beta) DNAJC5B CSPBETA
Observed MW:	22 kDa
Uniprot ID :	Q9UF47

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-DNAJC5B (IP:FNab02467, 3ug;
Detection:FNab02467 1:1500) with mouse testis
tissue lysate 4000ug.



mouse testis tissue were subjected to SDS PAGE
followed by western blot with
FNab02467(DNAJC5B antibody) at dilution of
1:1500