

DNAL4 antibody

Product Information

Catalog No.:	FNab02472
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Force generating protein of respiratory cilia. Produces force towards the minus ends of microtubules. Dynein has ATPase activity(By similarity).

Immunogen information

Immunogen:	dynein, axonemal, light chain 4
Synonyms:	Dynein axonemal light chain 4 DNAL4
Observed MW:	12 kDa
Uniprot ID :	O96015

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02472(DNAL4 antibody) at dilution of 1:300