

DNALI1 antibody

Product Information

Catalog No.:	FNab02473
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a dynamic role in flagellar motility.

Immunogen information

Immunogen:	dynein, axonemal, light intermediate chain 1
Synonyms:	Axonemal dynein light intermediate polypeptide 1 Inner dynein arm light chain, axonemal hp28 DNALI1
Observed MW:	30-33 kDa
Uniprot ID :	O14645

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis were subjected to SDS PAGE followed by western blot with FNab02473(DNALI1 antibody) at dilution of 1:5000