

DNASE1L1 antibody

Product Information

Catalog No.:	FNab02475
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	deoxyribonuclease I-like 1
Synonyms:	Deoxyribonuclease-1-like 1 DNase X Deoxyribonuclease I-like 1 (DNase I-like 1) Muscle-specific DNase I-like XIB DNASE1L1 DNAS1L1 DNL1L
Observed MW:	34 kDa
Uniprot ID :	P49184

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab02475(DNASE1L1 antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02475(DNASE1L1 antibody) at dilution of 1:1000