

DNASE1L3 antibody

Product Information

Catalog No.:	FNab02476
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Has DNA hydrolytic activity. Does not bind to actin. Cleaves chromatin DNA to nucleosomal units.

Immunogen information

Immunogen:	deoxyribonuclease I-like 3
Synonyms:	Deoxyribonuclease gamma (DNase gamma) DNase I homolog protein DHP2 Deoxyribonuclease I-like 3 (DNase I-like 3) Liver and spleen DNase (LS-DNase, LSD) DNASE1L3 DHP2 DNAS1L3
Observed MW:	33 kDa
Uniprot ID :	Q13609

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02476(DNASE1L3 antibody) at dilution of 1:1000