

TRDMT1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab02483 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|---|
| Immunogen: | tRNA aspartic acid methyltransferase 1 |
| Synonyms: | tRNA (cytosine(38)-C(5))-methyltransferase DNA (cytosine-5)-methyltransferase-like protein 2 (Dnmt2) DNA methyltransferase homolog HsaIIP (DNA MTase homolog HsaIIP, M.HsaIIP) PuMet TRDMT1 DNMT2 |
| Observed MW: | 43-45 kDa |
| Uniprot ID : | O14717 |

Application

| | |
|--------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image: | |



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02483(DNMT2 antibody) at dilution of 1:600