

DOCK1-Specific antibody

Product Information

Catalog No.:	FNab02493
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	dedicator of cytokinesis 1
Synonyms:	Dedicator of cytokinesis protein 1 180 kDa protein downstream of CRK (DOCK180) DOCK1
Observed MW:	180 kDa
Uniprot ID :	Q14185

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2400; IF: 1:10-1:100
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02493(DOCK180-Specific antibody) at dilution of 1:600