

## DXO antibody

### Product Information

Catalog No.:	FNab02507
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	dom-3 homolog Z(C. elegans)
Synonyms:	Decapping and exoribonuclease protein (DXO) 5'-3' exoribonuclease DXO Dom-3 homolog Z NAD-capped RNA hydrolase DXO (DeNADding enzyme DXO) DXO DOM3L DOM3Z NG6
Observed MW:	45 kDa
Uniprot ID :	O77932

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02507(DOM3Z antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02507(DOM3Z antibody) at dilution of 1:300