

FIZ1 antibody

Product Information

Catalog No.:	FNab02509
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	FLT3-interacting zinc finger 1
Synonyms:	Flt3-interacting zinc finger protein 1 Zinc finger protein 798 FIZ1 ZNF798
Observed MW:	60 kDa
Uniprot ID :	Q96SL8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab02509(domain-I-of-FIZ-1 Antibody) at dilution of 1:50



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab02509(domain-I-of-FIZ-1 antibody) at dilution of 1:200