

## DPCD antibody

### Product Information

Catalog No.:	FNab02512
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in the formation or function of ciliated cells.

### Immunogen information


Immunogen:	deleted in a Mouse model of primary ciliary dyskinesia
Synonyms:	Protein DPCD DPCD
Observed MW:	23-26 kDa
Uniprot ID :	Q9BVM2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart slide using FNab02512(DPCD Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab02512(RP11-529I10.4 antibody) at dilution of 1:300