

DPP7 antibody

Product Information

Catalog No.:	FNab02521
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays an important role in the degradation of some oligopeptides.

Immunogen information

Immunogen:	dipeptidyl-peptidase 7
Synonyms:	Dipeptidyl peptidase 2 Dipeptidyl aminopeptidase II Dipeptidyl peptidase 7 Dipeptidyl peptidase II (DPP II) Quiescent cell proline dipeptidase DPP7 DPP2 QPP
Observed MW:	54 kDa
Uniprot ID :	Q9UHL4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02521(DPP7 antibody) at dilution of 1:300