

DPPA4 antibody

Product Information

Catalog No.:	FNab02525
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the maintenance of active epigenetic status of target genes. May inhibit differentiation of embryonic cells into a primitive ectoderm lineage.

Immunogen information

Immunogen:	developmental pluripotency associated 4
Synonyms:	Developmental pluripotency-associated protein 4 DPPA4
Observed MW:	34-40 kDa
Uniprot ID :	Q7L190

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab02525(DPPA4 Antibody) at dilution of 1:300