

## DRD1 antibody

### Product Information

Catalog No.:	FNab02532
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Dopamine receptor whose activity is mediated by G proteins which activate adenylyl cyclase.

### Immunogen information

Immunogen:	dopamine receptor D1
Synonyms:	D(1A) dopamine receptor Dopamine D1 receptor DRD1
Observed MW:	50 kDa
Uniprot ID :	P21728

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02532(DRD1 antibody) at dilution of 1:300