

## DNM1L antibody

### Product Information

Catalog No.:	FNab02541
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	dynamin 1-like
Synonyms:	Dynamin-1-like protein Dnm1p/Vps1p-like protein (DVLP) Dynamin family member proline-rich carboxyl-terminal domain less (Dymple) Dynamin-like protein Dynamin-like protein 4 Dynamin-like protein IV (HdynIV) Dynamin-related protein 1 DNM1L DLP1 DRP1
Observed MW:	75-80 kDa
Uniprot ID :	O00429

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab02541(DRP1 antibody) at dilution of 1:100



Immunofluorescent analysis of HeLa cells using FNab02541 (DRP1 antibody) at dilution of 1:50



Various lysates were subjected to SDS PAGE followed by western blot with FNab02541(DRP1 antibody) at dilution of 1:1000