

DTL antibody

Product Information


Catalog No.:	FNab02549
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	denticleless homolog(Drosophila)
Synonyms:	Denticleless protein homolog DDB1- and CUL4-associated factor 2 Lethal(2) denticleless protein homolog Retinoic acid-regulated nuclear matrix-associated protein DTL CDT2 CDW1 DCAF2 L2DTL RAMP
Observed MW:	79 kDa
Uniprot ID :	Q9NZJ0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-DTL (IP:FNab02549, 3ug;
Detection:FNab02549 1:500) with mouse testis
tissue lysate 8000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab02549(DTL antibody) at
dilution of 1:200