

DTNA antibody

Product Information

Catalog No.:	FNab02550
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the formation and stability of synapses as well as being involved in the clustering of nicotinic acetylcholine receptors.

Immunogen information

Immunogen:	dystrobrevin, alpha
Synonyms:	Dystrobrevin alpha (DTN-A) Alpha-dystrobrevin Dystrophin-related protein 3 DTNA DRP3
Observed MW:	49-59 kDa, 78-85 kDa
Uniprot ID :	Q9Y4J8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab02550(DTNA Antibody) at dilution of 1:2000