

## DTWD2 antibody

### Product Information

Catalog No.:	FNab02552
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	DTW domain containing 2
Synonyms:	tRNA-uridine aminocarboxypropyltransferase 2 DTW domain-containing protein 2 DTWD2
Observed MW:	33 kDa
Uniprot ID :	Q8NBA8

### Application

Reactivity:	Human
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:2000
Image:	



IP Result of anti-DTWD2 (IP:FNab02552, 3ug;  
Detection:FNab02552 1:300) with COLO 320 cells  
lysate 2000ug.