

DUS4L antibody

Product Information

Catalog No.:	FNab02561
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the synthesis of dihydrouridine, a modified base found in the D-loop of most tRNAs.

Immunogen information

Immunogen:	dihydrouridine synthase 4-like
Synonyms:	tRNA-dihydrouridine(20a/20b) synthase [NAD(P)+]-like pp35 tRNA-dihydrouridine synthase 4-like DUS4L
Observed MW:	22, 30 kDa
Uniprot ID :	O95620

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02561(DUS4L antibody) at dilution of 1:1000