

DVL3 antibody

Product Information

Catalog No.:	FNab02574
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the signal transduction pathway mediated by multiple Wnt genes.

Immunogen information

Immunogen:	dishevelled, dsh homolog 3
Synonyms:	Segment polarity protein dishevelled homolog DVL-3 (Dishevelled-3) DSH homolog 3 DVL3 KIAA0208
Observed MW:	78 kDa
Uniprot ID :	Q92997

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab02574(DVL3 antibody) at dilution of 1:6000