

DCTN2 antibody

Product Information


Catalog No.:	FNab02581
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1B10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	dynactin 2(p50)
Synonyms:	Dynactin subunit 2 50 kDa dynein-associated polypeptide Dynactin complex 50 kDa subunit (DCTN-50) p50 dynamitin DCTN2 DCTN50
Observed MW:	50 kDa
Uniprot ID :	Q13561

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas slide using FNab02581(dynactin-2 Antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02581(dynactin-2 antibody) at dilution of 1:500