

DYNC2LI1 antibody

Product Information

Catalog No.:	FNab02587
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a motor for intraflagellar retrograde transport. Functions in cilia biogenesis(By similarity).

Immunogen information

Immunogen:	dynein, cytoplasmic 2, light intermediate chain 1
Synonyms:	Cytoplasmic dynein 2 light intermediate chain 1 (Dynein 2 light intermediate chain) DYNC2LI1 D2LIC LIC3
Observed MW:	40 kDa
Uniprot ID :	Q8TCX1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab02587(DYNC2LI1 antibody) at dilution of 1:500