

DYRK4 antibody

Product Information

Catalog No.:	FNab02593
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Possible non-essential role in spermiogenesis.

Immunogen information


Immunogen:	dual-specificity tyrosine-(Y)-phosphorylation regulated kinase 4
Synonyms:	Dual specificity tyrosine-phosphorylation-regulated kinase 4 DYRK4
Observed MW:	60 kDa
Uniprot ID :	Q9NR20

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab02593(DYRK4 antibody) at dilution of 1:100



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab02593(DYRK4 antibody) at dilution of 1:500