

SLC1A6 antibody

Product Information

Catalog No.:	FNab02607
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	solute carrier family 1(high affinity aspartate/glutamate transporter), member 6
Synonyms:	Excitatory amino acid transporter 4 Sodium-dependent glutamate/aspartate transporter Solute carrier family 1 member 6 SLC1A6 EAAT4
Observed MW:	60-70 kDa
Uniprot ID :	P48664

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab02607(EAAT4 antibody) at dilution of 1:300