

CDH1 antibody

Product Information

Catalog No.:	FNab02617
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cadherin 1, type 1, E-cadherin(epithelial)
Synonyms:	Cadherin-1 CAM 120/80 Epithelial cadherin (E-cadherin) Uvomorulin E-Cad/CTF1 E-Cad/CTF2 E-Cad/CTF3 CDH1 CDHE UVO
Observed MW:	125 kDa
Uniprot ID :	P12830

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC, FC
Recommended dilution:	WB: 1:5000-1:20000; IHC: 1:500-1:1000; IF: 1:20-1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab02617 (E-cadherin Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



Immunofluorescent analysis of (10% Formaldehyde) fixed A431 cells using FNab02617 (E-cadherin antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG (H+L)



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02617 (E-cadherin antibody) at dilution of 1:2000