

VOPP1 antibody

Product Information

Catalog No.:	FNab02628
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	EGFR-coamplified and overexpressed protein
Synonyms:	WW domain binding protein VOPP1 EGFR-coamplified and overexpressed protein (ECop) Glioblastoma-amplified secreted protein Putative NF-kappa-B-activating protein 055N Vesicular, overexpressed in cancer, prosurvival protein 1 VOPP1 ECOP GASP
Observed MW:	None
Uniprot ID :	Q96AW1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab02628(ECOP Antibody) at dilution of 1:50