

## EDC3 antibody

### Product Information

Catalog No.:	FNab02630
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds single-stranded RNA. Involved in the process of mRNA degradation and in the positive regulation of mRNA decapping. May play a role in spermiogenesis and oogenesis.

### Immunogen information


Immunogen:	enhancer of mRNA decapping 3 homolog
Synonyms:	Enhancer of mRNA-decapping protein 3 LSM16 homolog YjeF N-terminal domain-containing protein 2 (YjeF_N2, hYjeF_N2) YjeF domain-containing protein 1 EDC3 LSM16 YJDC YJEFN2
Observed MW:	56 kDa
Uniprot ID :	Q96F86

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab02630(EDC3 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02630(EDC3 antibody) at dilution of 1:500