

## ANP32B antibody

### Product Information

Catalog No.:	FNab02634
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	lysophosphatidic acid receptor 1
Synonyms:	Acidic leucine-rich nuclear phosphoprotein 32 family member B Acidic protein rich in leucines Putative HLA-DR-associated protein I-2 (PHAPI2) Silver-stainable protein SSP29 ANP32B APRIL PHAPI2
Observed MW:	45 kDa
Uniprot ID :	Q92633

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab02634(EDG2 Antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02634(EDG2 antibody) at dilution of 1:1000