

EDRF1 antibody

Product Information

Catalog No.:	FNab02638
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor involved in erythroid differentiation. Involved in transcriptional activation of the globin gene.

Immunogen information


Immunogen:	chromosome 10 open reading frame 137
Synonyms:	Erythroid differentiation-related factor 1 EDRF1 C10orf137
Observed MW:	140-160 kDa
Uniprot ID :	Q3B7T1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab02638(EDRF1 antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab02638(EDRF1 antibody) at dilution of 1:1000