

EEF1E1 antibody

Product Information

Catalog No.:	FNab02649
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Positive modulator of ATM response to DNA damage.

Immunogen information


Immunogen:	eukaryotic translation elongation factor 1 epsilon 1
Synonyms:	Eukaryotic translation elongation factor 1 epsilon-1 Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3 Elongation factor p18 Multisynthase complex auxiliary component p18 EEF1E1 AIMP3 P18
Observed MW:	20 kDa
Uniprot ID :	O43324

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02649(EEF1E1 antibody) at dilution of 1:10



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02649(EEF1E1 Antibody) at dilution of 1:1000