

## CRACR2A antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab02655  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | EF-hand calcium binding domain 4B   |
| Synonyms:    | EF-hand calcium-binding domain-containing protein 4B Calcium release-activated calcium channel regulator 2A (CRAC channel regulator 2A) Calcium release-activated channel regulator 2A Ras-related protein Rab-46 CRACR2A EFCAB4B RAB46 |
| Observed MW: | 46, 80-90 kDa   |
| Uniprot ID : | Q9BSW2  |

### Application

|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human        |
| Tested Application:      | ELISA, WB    |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image:                   |              |



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02655( EFCAB4B antibody at dilution of 1:600