

EFEMP2 antibody

Product Information

Catalog No.:	FNab02657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	EGF-containing fibulin-like extracellular matrix protein 2
Synonyms:	EGF-containing fibulin-like extracellular matrix protein 2 Fibulin-4 (FIBL-4) Protein UPH1 EFEMP2 FBLN4
Observed MW:	52-55 kDa
Uniprot ID :	O95967

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab02657(EFEMP2 antibody) at dilution of 1:50



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02657(EFEMP2 antibody) at dilution of 1:500