

EGR2 antibody

Product Information

Catalog No.:	FNab02672
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sequence-specific DNA-binding transcription factor. Binds to two specific DNA sites located in the promoter region of HOXA4. E3 SUMO-protein ligase helping SUMO1 conjugation to its coregulators NAB1 and NAB2, whose sumoylation down-regulates EGR2 own transcriptional activity.

Immunogen information

Immunogen:	early growth response 2(Krox-20 homolog, Drosophila)
Synonyms:	E3 SUMO-protein ligase EGR2 AT591 E3 SUMO-protein transferase ERG2 Early growth response protein 2 (EGR-2) Zinc finger protein Krox-20 EGR2 KROX20
Observed MW:	50 kDa
Uniprot ID :	P11161

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02672(EGR2 antibody) at dilution of 1:1000