

EI24 antibody

Product Information

Catalog No.:	FNab02680
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a negative growth regulator via p53-mediated apoptosis pathway. Regulates formation of degradative autolysosomes during autophagy(By similarity).

Immunogen information

Immunogen:	etoposide induced 2.4 mRNA
Synonyms:	Etoposide-induced protein 2.4 homolog p53-induced gene 8 protein EI24 PIG8
Observed MW:	40-45 kDa
Uniprot ID :	O14681

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02680(EI24 antibody) at dilution of 1:600