

EIF1 antibody

Product Information

Catalog No.:	FNab02683
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Necessary for scanning and involved in initiation site selection. Promotes the assembly of 48S ribosomal complexes at the authentic initiation codon of a conventional capped mRNA.

Immunogen information


Immunogen:	eukaryotic translation initiation factor 1
Synonyms:	Eukaryotic translation initiation factor 1 (eIF1) A121 Protein translation factor SUI1 homolog Sui1 iso1 EIF1 SUI1
Observed MW:	13 kDa
Uniprot ID :	P41567

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using FNab02683(EIF1 Antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab02683(EIF1 Antibody) at dilution of 1:1000