

EIF1B antibody

Product Information

Catalog No.:	FNab02688
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probably involved in translation.

Immunogen information


Immunogen:	eukaryotic translation initiation factor 1B
Synonyms:	Eukaryotic translation initiation factor 1b (eIF1b) Protein translation factor SUI1 homolog GC20 EIF1B
Observed MW:	14 kDa
Uniprot ID :	O60739

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab02688(EIF1B Antibody) at dilution of 1:50



L02 cells were subjected to SDS PAGE followed by western blot with FNab02688(EIF1B antibody) at dilution of 1:300