

EIF2B1 antibody

Product Information

Catalog No.:	FNab02693
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

Immunogen information

Immunogen:	eukaryotic translation initiation factor 2B, subunit 1 alpha, 26kDa
Synonyms:	Translation initiation factor eIF2B subunit alpha eIF2B GDP-GTP exchange factor subunit alpha EIF2B1 EIF2BA
Observed MW:	34 kDa
Uniprot ID :	Q14232

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02693(EIF2B1 antibody) at dilution of 1:500