

EIF2B2 antibody

Product Information

Catalog No.:	FNab02694
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

Immunogen information


Immunogen:	eukaryotic translation initiation factor 2B, subunit 2 beta, 39kDa
Synonyms:	Translation initiation factor eIF2B subunit beta S20I15 S20III15 eIF2B GDP-GTP exchange factor subunit beta EIF2B2 EIF2BB
Observed MW:	39 kDa
Uniprot ID :	P49770

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab02694(EIF2B2 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02694(EIF2B2 Antibody) at dilution of 1:1000