

## EIF2B4 antibody

### Product Information

Catalog No.:	FNab02696
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

### Immunogen information


Immunogen:	eukaryotic translation initiation factor 2B, subunit 4 delta, 67kDa
Synonyms:	Translation initiation factor eIF2B subunit delta eIF2B GDP-GTP exchange factor subunit delta EIF2B4 EIF2BD
Observed MW:	60 kDa
Uniprot ID :	Q9UI10

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02696(EIF2B4 antibody) at dilution of 1:100



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02696(EIF2B4 antibody) at dilution of 1:400