

AGO3-Specific antibody

Product Information

Catalog No.:	FNab02698
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	eukaryotic translation initiation factor 2C, 3
Synonyms:	Protein argonaute-3 (Argonaute3, hAgo3) Argonaute RISC catalytic component 3 Eukaryotic translation initiation factor 2C 3 (eIF-2C 3, eIF2C 3) AGO3 EIF2C3
Observed MW:	95 kDa
Uniprot ID :	Q9H9G7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab02698(EIF2C3-Specific Antibody) at dilution of 1:300