

EIF3H antibody

Product Information

Catalog No.:	FNab02709
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	eukaryotic translation initiation factor 3, subunit H
Synonyms:	Eukaryotic translation initiation factor 3 subunit H (eIF3h) Eukaryotic translation initiation factor 3 subunit 3 eIF-3-gamma eIF3 p40 subunit EIF3H EIF3S3
Observed MW:	40 kDa
Uniprot ID :	O15372

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200; IP: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNa02709(EIF3H antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNa02709(EIF3H antibody) at dilution of 1:1000