

## EIF4E3 antibody

### Product Information

Catalog No.:	FNab02722
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Recognizes and binds the 7-methylguanosine-containing mRNA cap during an early step in the initiation of protein synthesis. May act as an inhibitor of EIF4E1 activity(By similarity).

### Immunogen information

Immunogen:	eukaryotic translation initiation factor 4E family member 3
Synonyms:	Eukaryotic translation initiation factor 4E type 3 (eIF-4E type 3, eIF-4E3, eIF4E type 3, eIF4E-3) EIF4E3
Observed MW:	24 kDa
Uniprot ID :	Q8N5X7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02722(EIF4E3 antibody) at dilution of 1:500