

EMILIN1 antibody

Product Information

Catalog No.:	FNab02755
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be responsible for anchoring smooth muscle cells to elastic fibers, and may be involved not only in the formation of the elastic fiber, but also in the processes that regulate vessel assembly. Has cell adhesive capacity.

Immunogen information


Immunogen:	elastin microfibril interfacer 1
Synonyms:	EMILIN-1 Elastin microfibril interface-located protein 1 (Elastin microfibril interfacer 1) EMILIN1 EMI
Observed MW:	120 kDa
Uniprot ID :	Q9Y6C2

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab02755(EMILIN1 antibody) at dilution of 1:100



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab02755(EMILIN1 antibody) at dilution of 1:300