

ENPP5 antibody

Product Information

Catalog No.:	FNab02775
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in neuronal cell communication. Lacks nucleotide pyrophosphatase and lysopholipase D activity(By similarity).

Immunogen information


Immunogen:	ectonucleotide pyrophosphatase/phosphodiesterase 5(putative function)
Synonyms:	Ectonucleotide pyrophosphatase/phosphodiesterase family member 5 (E-NPP 5, NPP-5) ENPP5
Observed MW:	55 kDa
Uniprot ID :	Q9UJA9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab02775(ENPP5 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02775(ENPP5 antibody) at dilution of 1:300