

ENTPD7 antibody

Product Information

Catalog No.:	FNab02786
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Preferentially hydrolyzes nucleoside 5'-triphosphates. The order of activity with respect to possible substrates is UTP > GTP > CTP.

Immunogen information

Immunogen:	ectonucleoside triphosphate diphosphohydrolase 7
Synonyms:	Ectonucleoside triphosphate diphosphohydrolase 7 (NTPDase 7) Lysosomal apyrase-like protein 1 ENTPD7 LALP1
Observed MW:	69 kDa
Uniprot ID :	Q9NQZ7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02786(ENTPD7 Antibody) at dilution of 1:300