

EVPL antibody

Product Information

Catalog No.:	FNab02789
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	envoplakin
Synonyms:	Envoplakin 210 kDa cornified envelope precursor protein 210 kDa paraneoplastic pemphigus antigen p210 EVPL
Observed MW:	210 kDa
Uniprot ID :	Q92817

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab02789(Envoplakin antibody at dilution of 1:300