

RAPGEF4 antibody

Product Information

Catalog No.:	FNab02793
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Rap guanine nucleotide exchange factor(GEF) 4
Synonyms:	Rap guanine nucleotide exchange factor 4 Exchange factor directly activated by cAMP 2 Exchange protein directly activated by cAMP 2 (EPAC 2) cAMP-regulated guanine nucleotide exchange factor II (cAMP-GEFII) RAPGEF4 CGEF2 EPAC2
Observed MW:	110 kDa
Uniprot ID :	Q8WZA2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02793(EPAC2 antibody) at dilution of 1:1000