

EPB41L2 antibody

Product Information

Catalog No.:	FNab02797
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	erythrocyte membrane protein band 4.1-like 2
Synonyms:	Band 4.1-like protein 2 Erythrocyte membrane protein band 4.1-like 2 Generally expressed protein 4.1 (4.1G) EPB41L2
Observed MW:	130-150 kDa
Uniprot ID :	O43491

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02797(EPB41L2 antibody) at dilution of 1:800