

## EPS15L1 antibody

### Product Information

Catalog No.:	FNab02818
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	epidermal growth factor receptor pathway substrate 15-like 1
Synonyms:	Epidermal growth factor receptor substrate 15-like 1 Eps15-related protein (Eps15R) EPS15L1 EPS15R
Observed MW:	94 kDa
Uniprot ID :	Q9UBC2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab02818(EPS15R antibody) at dilution of 1:500