

ERMAP antibody

Product Information

Catalog No.:	FNab02850
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Possible role as a cell-adhesion or receptor molecule of erythroid cells.

Immunogen information

Immunogen:	erythroblast membrane-associated protein(Scianna blood group)
Synonyms:	Erythroid membrane-associated protein (hERMAP) Radin blood group antigen Scianna blood group antigen ERMAP RD SC
Observed MW:	53 kDa
Uniprot ID :	Q96PL5

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-ERMAP (IP:FNab02850, 4ug;
Detection:FNab02850 1:300) with MCF-7 cells
lysate 1600ug.