

## ESAM antibody

### Product Information

Catalog No.:	FNab02860
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Can mediate aggregation most likely through a homophilic molecular interaction.

### Immunogen information

Immunogen:	endothelial cell adhesion molecule
Synonyms:	Endothelial cell-selective adhesion molecule ESAM
Observed MW:	56-58 kDa
Uniprot ID :	Q96AP7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human lung tissue were subjected to SDS PAGE followed by western blot with FNab02860(ESAM antibody) at dilution of 1:500